

This application is in accordance with [Regulation no. 200/2020](#) on the Importation of Dogs and Cats. Information on importation (health conditions and quarantine) are on the [MAST website](#) and it is recommended to read the instructions before applying for an import permit. Dogs can only be imported to Iceland from approved exporting countries¹⁾

Applicant / Importer

Name	Tel.no.
Address	
e-mail:	Country

Dog

Dog breed ²⁾	Country of export
Date of birth	Name of dog
Microchip number ³⁾	

By submitting this application the applicant confirms to comply in every respect with the conditions of the importation of a dog according to regulation no. 200/2020 on the Importation of Dogs and Cats. Dogs that do not fully comply with import conditions will be exported or euthanised at the expense of the importer.

Place

Date

Signature of applicant

A scanned copy of the signed application must be sent by e-mail to petimport@mast.is

Applicants with Icelandic ssn (kennitala) should apply on the [MAST portal](#).

1) Approved exporting countries

Category 1 countries (rabies free countries):

Australia, Belgium, Britain, Bulgaria, Croatia, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Ireland, Italy, Japan, Latvia, Liechtenstein, Luxembourg, Malta, Netherlands, Norway (excl. Svalbard), Portugal, Slovakia, Slovenia, Spain, Sweden, Switzerland, The Faroe Islands

Category 2 countries (rabies well controlled):

Canada, Greenland, Lithuania, Poland, Romania, Serbia, Singapore, Turkey, Thailand, USA

2) Prohibited dog breeds.

The following breeds are not permitted in Iceland:

Hybrids from wolves and dogs for at least 10 generations. The following dog breeds, as well as their cross breeds for at least 5 generations. 1) Pit Bull Terrier/Staffordshire Bull Terrier, 2) Fila Brasileiro, 3) Tosa Inu, 4) Dogo Argentino. MAST can reject other breeds considered to be a risk. Dog breeds that are closely related to the above mentioned breeds have been rejected by MAST. If in doubt about the type of dog you are planning to import, please contact MAST for further information.

3) Microchip must be ISO compliant